

See how the Personal Medik8 Titanium Dermastamp can help with your skin problem:

Lip Wrinkles

Wrinkles around the mouth can become more pronounced due to facial expressions such as smiling, frowning and even smoking, promoting the formation of fine lines around the lips.

Piercing the upper layers of the skin with microneedles leads to increased blood flow in the skin and stimulates the synthesis of collagen and elastin. Enhanced penetration of topically applied Medik8 products will also assist in stimulating the production of these skin proteins. Each subsequent treatment gradually reduces the appearance of lines and wrinkles.

Use the 0.2mm Titanium Dermastamp if you have fine lines or very thin skin around the lips, or if you are just beginning your home-based Skin Needling treatment.

Use the 0.3mm Titanium Dermastamp if your skin around the lips is comparatively thick or if you have deeper peri-oral lines.

Notice: To build up skin tolerance to skin needling, allow 3 day intervals between treatments for a minimum of 2 weeks. As your skin becomes adjusted to needling, gradually increase the frequency of use. Do not use more than once a day. If used more than once a day the Medik8 Titanium Dermastamp may cause more inflammation than is aimed for and even cause damage.

Medik8 Titanium Dermastamp Lip Wrinkle Protocol - For the Lip Area

1. **STERILISE:** Disinfect the Titanium Dermastamp by using Medik8 Sanitising Solution directly on the needles before every use and allow to dry.
2. **CLEANSE:** Cleanse the whole peri-oral area concerned, either above the lips and in the naso-labial folds and/ or below the lips, by applying Medik8 betaCleanse or your preferred Medik8 cleanser, following product instructions. Ensure the skin is left clean and dry.
3. **STAMP:** Hold the Dermastamp (like a pen) close to the needle head, providing greater control and accuracy of application. The free hand can be used to hold the skin taught surrounding the lips to stretch the lines/ wrinkles.

Rest your elbow on a flat surface to avoid movement of the needles when in contact with your skin.

Hold the Dermastamp directly over the targeted wrinkle. Repeatedly press the needles into the skin 20-25 times ensuring the whole length of the wrinkle receives repetitive needling.

IMPORTANT TIP: fully remove the Dermastamp from contact with the skin between each application. Do not rotate/ twist the Dermastamp whilst in contact with the skin.

4. **PREVENT:** Apply Medik8 vitamin C serum: C-Tetra or CE-Tetra directly over the stamped area and allow to fully absorb.
5. **CORRECT:** Apply Medik8 Retinol 3 TR directly over the stamped area and allow to fully absorb.

Optional: Apply Pretox pout to the lips for a complete lip rejuvenation. Pretox pout is also effective when applied beyond the lip line, directly onto the stamped area, to optimise the lip anti-ageing protocol.*

6. FINISH: Clean the Dermastamp by running under hot water and sanitise as before.

Allow to dry before replacing cap.

*Application of these products should proceed the application of vitamin C and Retinol serums in line with the Medik8 protocol.

Treatment duration: 5 – 10 minutes

Frequency of use: 2 -7 times per week, depending on skin tolerance.

To build up skin tolerance to Skin Needling, allow 3 day intervals between treatments for a minimum of 2 weeks. Gradually increase frequency of use, when your skin becomes adjusted to needling.

Do not use more than once a day.

Optimal treatment time: For best results, use at night to allow for maximum absorption of the active formulation