

FAQ's - SONIC SPATULA

How can ultrasonic waves cleanse my skin?

The Sonic Spatula generates high frequency sonic vibration (ultrasonic waves) that cause the skin cells to vibrate and the tissue to become softened. The water on the skin is pushed with the ultrasonic waves creating micro-jets which extract environmental toxins, makeup and other impurities that are in the pores, thereby effectively cleansing and exfoliating the skin.

How can ultrasonic waves also penetrate products?

The speed of the vibration and the tip of the Sonic Spatula direct the ultrasonic waves. So when in contact with a water-based product like gels or serums, the ultrasonic wave pushes in the product for greater penetration and efficacy.

Does the SonicSpatula hurt or sting the skin?

The Ultrasonic vibration waves do not cause any uncomfortable sensation or hurt the skin. It is important to have either water or a product on the skin so the tip can glide.

Why should I use Rosewater or Lavender Misting Spray with the SonicSpatula rather than regular water?

By using the Energized Misting Sprays with the Sonic Spatula, the vibration and water create "micro-jets" which thoroughly exfoliate the skin. Made with magnetized and vortexed water, which is believed to raise the energy to a higher vibration, the misting sprays add a therapeutic benefit to your treatment.

Can I use anything other than TEI Spa Misting Sprays with my Sonic Spatula?

Yes, you may use any toner, cleanser or misting spray when using the exfoliation setting of the Sonic Spatula. The main concern is the ability for the spatula head to glide across the skin and cause no discomfort.

What skin types should use the SonicSpatula?

This tool is designed for all skin types. Dry skin can benefit greatly from using the Sonic Spatula by removing dead skin and aiding in the penetration of products. Normal to oily skin types can achieve brighter and clearer skin by using the Sonic Spatula on a regular basis. Combination skin can use the Sonic Spatula as needed for deep cleansing and product penetration. Please refer to the table for additional information.

Skin Type - Oily Skin - Duration of Use - 3-4 min – Frequency - Up to 3-4 times per week

Skin Type - Combination Skin - Duration of Use - 2-3 min – Frequency - Up to 2-3 times per week

Skin Type - Dry/Sensitive Skin - Duration of Use -1-2 min – Frequency - Up to 1-2 times per week